

Product Datasheet

Mammal IL-31 protein (orb1215822)

Description	The Dolphin IL-31 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	IL-31
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SHTAPYVLGP DKQQAIIKEL RSSFEKLLND YQQEESGVPE SSQYQMPCFA LNPQLPYNID SSATLPYCRV IKPLFDNKAQ VVALIDHLSK LEFRHEPETK VSVPTDSFKC KSFIVSIFKK FSECLICLNH (130)
Purity	98%
Source	Yeast
MW	14.8 kDa
Entrez	101334968
Expiration Date	6 months from date of receipt.