

Product Datasheet

Rat CD40LG protein (orb1215831)

Description	The Rat CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQKGDEDPQI AAHVVSEANS NAASVLQWAK KGYTMSNL VVLENGRQLT VKREGLYYVY TQVTFCSNRE PLSQRPFIVS LWLKPSSGSE RILLRAANTH SSSKLCEQQS IHLGGVFELQ AGASVFNVT EASQVIHGIG FSSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.4 kDa
Entrez	84349
Expiration Date	6 months from date of receipt.