

Product Datasheet

Monkey IFN alpha protein (orb1215836)

| | |
|-------------------------|---|
| Description | The Cynomolgus Monkey IFN alpha 2 yeast-derived recombinant protein is not tagged and is naturally e... |
| Target | IFN alpha |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | CDLPQTHSLG NRRTLMLLAQ MRRISLFFCL KDRHDFEFPQ EEFGNQFQKA QTIPVLHEMI QQTFNLFSTK DSSAAWDETL LNKFYTELYQ QLNDLEACVM QEMGVTETPL MNKNSILAVR KYFQRITLYL KKKYSLCAW EVVRAEIMRS FSLSTNLQES LRSKE (165) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.5 kDa |
| Entrez | 102142148 |
| Expiration Date | 6 months from date of receipt. |