

## Product Datasheet

### Porcine Midkine protein (orb1215841)

<b>Description</b>	The Swine Midkine yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	Midkine
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	KKKDKVKKGG PGSECAEWTW GPCTPSSKDC GVGFREGTCG AQTQRVRCRV PCNWKKEFGA DCKYKFESWG ACDGGTGTKA RQGLKKARY NAQCQETIRV TKPCTPKTKA KAKAKKGK GK D (121)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	13.2 kDa
<b>Entrez</b>	<b>100359358</b>
<b>Expiration Date</b>	6 months from date of receipt.