

Product Datasheet

Rabbit Midkine protein (orb1215845)

| | |
|-------------------------|---|
| Description | The Rabbit Midkine yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ... |
| Target | Midkine |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | KKKDKVKGGP GSECAEWTWG PCTPSSKDCG VGFREGTCGA QTQRIRCRLP CNWKKEFGAD CKYKFDSWGA CDGGTGTKAR QGTLKKARYN AQCQETIRVT KPCTPKTKAK AKAKKGKAKD (120) |
| Purity | 98% |
| Source | Yeast |
| MW | 13.1 kDa |
| Entrez | 100008929 |
| Expiration Date | 6 months from date of receipt. |