

Product Datasheet

Ferret IL-1 alpha protein (orb1215850)

Description	The Ferret IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SAAYNFRSNE KYNMRIKS QFILNDNLNQ SIIQQTGGNY LMAAALQNL DAVKFDMGAY TSEDSKFPVT LRISKTRLFV SAQNEGEPVL LKEMPETPKT IKDETLLFF WERHGTKNYF TSVAQPKLFI ATQERKLVHM AKGQPSITDF QILET (155)
Purity	98%
Source	Yeast
MW	17.8 kDa
Entrez	101691850
Expiration Date	6 months from date of receipt.