

## Product Datasheet

### Rabbit Amphiregulin protein (orb1215865)

<b>Description</b>	The Rabbit Amphiregulin yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
<b>Target</b>	Amphiregulin
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SVRVEQVIKP KKNKAESEKT SDKPKRKKKG GKNGKNRKNK KKKNPCDAEF QNFCIHGECK YIEHLDAVTC KCHHDYFGER CGEKSMK (87)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	10.1 kDa
<b>Entrez</b>	<a href="#">100354550</a>
<b>Expiration Date</b>	6 months from date of receipt.