

## Product Datasheet

### Human Prolactin protein (orb1215877)

|                         |   |
|-------------------------|---|
| <b>Description</b>      | The Human Prolactin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....   |
| <b>Target</b>           | Prolactin   |
| <b>Form/Appearance</b>  | Lyophilized   |
| <b>Storage</b>          | -20°C   |
| <b>Note</b>             | For research use only   |
| <b>Protein Sequence</b> | LPICPGGAAR CQVTLRDLFD RAVVLSHYIH<br>NLSSEMFSEF DKRYTHGRGF ITKAINSCHT<br>SSLATPEDKE QAQQMNQKDF LSLIVSILRS<br>WNEPLYHLVT EVRGMQEAPE AILSKAVEIE<br>EQTKRLLGEM ELIVSQVHPE TKENEIYPVW<br>SGLPSLQ MAD EESRLSAYYN LLHCLRRDSH<br>KIDNYLKLLK CRIIHNNNC (199) |
| <b>Purity</b>           | 98%   |
| <b>Source</b>           | Yeast   |
| <b>MW</b>               | 22.9 kDa  |
| <b>Entrez</b>           | <b>5617</b>   |
| <b>Expiration Date</b>  | 6 months from date of receipt.  |