

Product Datasheet

Zebrafish CCL3 protein (orb1215883)

Description	The Zebrafish CCL3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	CCL3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	DVMGNHANTP DACCFTEFKR KIHPSKINSY NPTRVDCTLP GVIFVTQKGL RLCVEPKLNW VKKTIQIIDD RNI (73)
Purity	98%
Source	Yeast
MW	8.4 kDa
Entrez	100537088
Expiration Date	6 months from date of receipt.