

Product Datasheet

Porcine OPN protein (orb1215886)

Description	The Swine Osteopontin yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	OPN
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	LPVKQTNSGS SEEKLLSNKY TDAVATLLKP DPSQKQTFLA PQNTISSEET DDFKQETLPSKSNESPEQTD DVDDDDDEDH VDSRDTDSEE ADHADDADRS DESHHSDESD ELVTDFPTDTPATDVTPAVP TGDPNDRGRGD SVVYGLRSKS KKFRRSEAQQ LDATEEDLTS HVESEETDGTPKAILVAQRL HVASDLDSQE KDSQETSQPD DRSVETRSE QSKEYTIKTY DGSNEHSNVIESQENPKVSQ EFHSHEDKLV PDSKSEEDKH LKLRSHELE SASSEIN (287)
Purity	98%
Source	Yeast
MW	31.9 kDa
Expiration Date	6 months from date of receipt.