

Product Datasheet

Porcine CTLA-4 protein (orb1215889)

Description	The Swine sCTLA-4 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	CTLA-4
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MHVAQPAVVL ANSRGVASFV CEYGSAGKAA EVRVTVLRRA GSQMTEVCAA TYTVEDELTF LDDSTCTGTS TENKVNLTIQ GLRAVDTGLY ICKVELLYPP PYYVGMNGT QIYVIDPEPC PDS (124)
Purity	98%
Source	Yeast
MW	13.3 kDa
Entrez	397286
Expiration Date	6 months from date of receipt.