

Product Datasheet

Human CX3CL1 protein (orb1215896)

Description	The Human CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th...
Target	CX3CL1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QHHGVTKCNI TCSKMTSKIP VALLIHYQQN QASCGKRAII LETRQHRLFC ADPKEQWVKD AMQHLDQRQAA ALTRNG (76)
Purity	98%
Source	Yeast
MW	8.6 kDa
Entrez	6376
Expiration Date	6 months from date of receipt.