

Product Datasheet

Mouse TIM-3 protein (orb1215909)

Description	The Mouse TIM-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	TIM-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	RSLENAYVFE VGKNAYLPCS YTLSTPGALV PMCWKGKFCP WSQCTNELLR TDERNVTYQK SSRYQLKGD L NKGDVSLIK NVTLDDHGTY CCRIQFPGLM NDKKLELKD IKAQVTPAQ TAHGDSTTAS PRTLTERNG SETQTLVTLH NNNGTKISTW ADEIKDSGET IRTA (174)
Purity	98%
Source	Yeast
MW	19.3 kDa
Entrez	171285
Expiration Date	6 months from date of receipt.