

Product Datasheet

Mouse Ang-1 protein (orb1215910)

Description	The Mouse Angiogenin-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	Ang-1
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	QDDSRVTKFL TQHHDAPKPG RDDRYCERMM KRRSLTSPCK DVNTFIHGK SNIKAICGAN GSPYRENLRM SKSPFQVTTTC KHTGGSPRPP CQYRASAGFR HVVIACENGL PVHFDESFFS L (121)
Purity	98%
Source	Yeast
MW	13.7 kDa
Entrez	11727
Expiration Date	6 months from date of receipt.