

## Product Datasheet

### Rabbit Ang-1 protein (orb1215918)

<b>Description</b>	The Rabbit Angiogenin-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-f...
<b>Target</b>	Ang-1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	QDDSRYPKHFL TQHYDAKPFG RNDRYCETMM KRRDLTSPCK DTNTFVHGK GSIKDVCEK NGKPYGKNFR ISKSSFQVTT CKHVGGSPWP PCRYRATSGS RNIVACENG LPVHFDESVF QKKAH (125)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.3 kDa
<b>Entrez</b>	<b>100354306</b>
<b>Expiration Date</b>	6 months from date of receipt.