

Product Datasheet

Human EGF protein (orb1215922)

Description	The Human EGF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	EGF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	NSDSECPLSH DGYCLHDGVC MYIEALDKYA CNCVVG YIGE RCQYRDLKWW ELR (53)
Purity	98%
Source	Yeast
MW	6.2 kDa
Entrez	1950
Expiration Date	6 months from date of receipt.