

Product Datasheet

Fish IL-4-13 protein (orb1215929)

Description	The Atlantic Salmon IL-4-13 yeast-derived recombinant protein is not tagged and is naturally endotox...
Target	IL-4-13
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	YRQHKNHHSN LEAIISMAEQ YTHNLTEDLL RTLLIDVTHL TGTTTKCKEF FCEAEKILAS VKNATFGNIV RHLRVYNKHQ YCIVVKSQDN QAEGNQVELR RLLLDLKECG QKINSKP (117)
Purity	98%
Source	Yeast
MW	13.5 kDa
Entrez	100534607
Expiration Date	6 months from date of receipt.