

Product Datasheet

Human PDGF-BB protein (orb1215937)

Description	The Human PDGF-BB yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
Target	PDGF-BB
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SLGSLTIAEP AMIAECKTRT EVFEISRRLI DRTNANFLVW PPCVEVQRCS GCCNNRNVQC RPTQVQLRPV QVRKIEIVRK KPIFKKATVT LEDHLACKCE TVAAARPVT (109)
Purity	98%
Source	Yeast
MW	12.3 kDa (monomer)
Entrez	5155
Expiration Date	6 months from date of receipt.