

Product Datasheet

Fish IL-6 protein (orb1215949)

Description	The Atlantic Salmon IL-6 yeast-derived recombinant protein is not tagged and is naturally endotoxin-...
Target	IL-6
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SPVPSALSEL MTSGWTSGEE LATDSLYGET GAPPKWEQII KMLVHEVTTL RNQQFVEEFQ KPVVEISSFS QHQLVSTPPH LSKTLCSTTN KEACLQEIRR GLRVYQLLLQ HVKAEYPOST LLPSVTHQTT VLIQLVKDQM KAAEVVEDLS ASERERVLGE VSTGTEWERK TSVHAILREL RNFLVDTKRA LRRMEKRGKD FQ (202)
Purity	98%
Source	Yeast
MW	22.9 kDa
Entrez	106570581
Expiration Date	6 months from date of receipt.