

Product Datasheet

Porcine IL-3 protein (orb1215972)

Description	The Swine IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-3
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MPTTTLQPKN YLAMIQEITR SLENLTVTSN KSLTLNELET LVNNTLLRPN LEAFVTF AEN HLKNISGIKK NLEKFRPILP TSMSTEEPIS IEEGDLGDFR AKLMEYLVVL RDSLKPMITE P (121)
Purity	98%
Source	Yeast
MW	13.7 kDa
Expiration Date	6 months from date of receipt.