

Product Datasheet

Mouse CD40LG protein (orb1215978)

Description	The Mouse CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQRGDEDPQI AAHVVSEANS NAASVLQWAK KGYTMSKNL VMLENGKQLT VKREGLYYVY TQVTFCSNRE PSSQRPFIG LWLKPSSGSE RILLKAANTH SSSQLCEQQS VHLGGVFELQ AGASVFNVT EASQVIHRVG FSSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.4 kDa
Entrez	21947
Expiration Date	6 months from date of receipt.