

## Product Datasheet

### Bovine IL-3 protein (orb1215982)

<b>Description</b>	The Bovine IL-3 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	IL-3
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	APQAKGLPVV TS RTPYSMLM KEIMDDLKKI TPSPEGLNS DEKNFLTKES LLQANLKVFM TFATDTFGSD SKIMKNLKEF QPVLPTATPT EDPIFIENKN LGDFRMKLEE YLVIIRNYLK SKNIWFS (127)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.5 kDa
<b>Entrez</b>	<b>280823</b>
<b>Expiration Date</b>	6 months from date of receipt.