

Product Datasheet

Rabbit PD-L1 protein (orb1215984)

| | |
|-------------------------|---|
| Description | The Rabbit PD-L1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | PD-L1 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | FTVTVPKDLY VVEYGSNVTI ECKFPVEKQL DLNSLVVYWE MGDKQIIQFV HGNENLKVQH SSYRQRAQLL KDQLHLGIAA LQITDVKLRD AGVYCCMISY GGADYKRITL KVNAPYHKIN QRISVDPVTS EHELTCAEG YPEAEVIWTS SDHQILHGKT DVTNSNREEK FFNVTSTLRI NTTANEIFYC TFQRPGEAN NTAKLIPEP PVFHRNER (218) |
| Purity | 98% |
| Source | Yeast |
| MW | 24.9 kDa |
| Entrez | 100357106 |
| Expiration Date | 6 months from date of receipt. |