

Product Datasheet

Porcine CBG protein (orb1215995)

Description	The Swine CBG yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The p...
Target	CBG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	EENFYVNETA TVKVPMMFQS RAMKYLNDL LPCQLVQLEY TGNETAFFIL PVKGEMDTVI AGLSRDTIQR WSKSLIPSQV DLYVPKVSIS GAYDLGSILG DMGIVDLLSH PTHFSGITQN ALPKMSKVH KAVLQFDEKG MEAAAPTTRG RSLHAAPKPV TVH (163)
Purity	98%
Source	Yeast
MW	17.9 kDa
Entrez	396736
Expiration Date	6 months from date of receipt.