

## Product Datasheet

### Monkey IL-4 protein (orb1216016)

<b>Description</b>	The Cynomolgus Monkey IL-4 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
<b>Target</b>	IL-4
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	HKCDITLQEI IKTLNSLTEQ KTLCTKLTIT DILAASKNTT EKETFCRAAT VLRQFYSHHE KDTRCLGATA QQFHRHKQLI RFLKRLDRNL WGLAGLNSCP VKEANQSTLE NFLERLKTIM REKYSKCSS (129)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	14.9 kDa
<b>Entrez</b>	<b>107130068</b>
<b>Expiration Date</b>	6 months from date of receipt.