

Product Datasheet

Porcine CXCL10 protein (orb1216019)

Description	The Cynomolgus Monkey CXCL10 yeast-derived recombinant protein is not tagged and is naturally endoto...
Target	CXCL10
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IPLSRTVRCT CISISNPVN PRSLEKLEII PPSQFCPHVE IIATMKKKGE KRCLNPESKA IKNLLKAVSK ERSKRSP (77)
Purity	98%
Source	Yeast
MW	8.7 kDa
Entrez	102134351
Expiration Date	6 months from date of receipt.