

Product Datasheet

Goat IFN alpha protein (orb1216023)

Description	The Caprine IFN alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IFN alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	ICSLGCHLPH IHSLANGRVL MLLRQLRRVS PSSCLQDRKD FLFPQEALGG SQLQKAQALS VLHEVTQHTF QLFSTEGTAA AWDQSLLDKL RAALDQQLTD LQACLRQEEG LRGAPLLKED SSLAVRKYFH RVTLYLQKEG HSPCAWEVVR AEVMRAFSSS TNLQERFRRK D (171)
Purity	98%
Source	Yeast
MW	19.4 kDa
Entrez	100860850
Expiration Date	6 months from date of receipt.