

Product Datasheet

Sheep IL-12 p40 protein (orb1216027)

Description	The Ovine/Caprine IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally e...
Target	IL-12 p40
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IWELEKNVYV VELDWYPNAP GETVVLTCDT PEEDGITWTS DQSSEVLGSG KTLTIQVKEF GDAGQYTCHK GGEVLSRSL LLHKKEDGIW STDILKDQKE PKAKSFLKCE AKDYSGHFTC SWLTAISTNL KFSVKSSRGS SDPRGVTCGA ASLSAEKVSM DHREYNKYTV ECQEGSACPA AEESLPIEVV MEAVHKLKYE NYTSSFFIRD IIKPDPPKNL QLRPLKNSRQ VEVSWEYPDT WSTPHSYFSL TFCVQVQGKN KREKKLFTDQ TSAKVTCHKD ANIRVQARDR YYSSFWSEWA SVSCS (305)
Purity	98%
Source	Yeast
MW	34.5 kDa
Entrez	443472
Expiration Date	6 months from date of receipt.