

## Product Datasheet

### Bovine ICAM-1 protein (orb1216037)

<b>Description</b>	The Bovine sICAM-1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
<b>Target</b>	ICAM-1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GISIHPSKAI IPRGDSLTVN CSNSCDQKST FGLETVLIKE EVGRGDNWKV FQLRDVQEDI ELFCYSNCHK EQTIASMNLT VYWFPEHVEL APLPLWQPVG EELNLSCLVS GGAPRAHLSV VLLRGEELG RQPVGKGEPA KVMFTVQSRR EDHGTNFSCR WELDLRSQGL ELFQNTSAPR KLQTYVLPSI DPHLEVPIV EVGSRWPVNC TLDGLFPASD AKVYLVLDQ KLESNITYDG DSVLAKAWME ENEEGTHSLK CSVTLGEVSR RTQENVTVYS FPLPTLTLSP PEVSEWTTVT VECVTRDGAV VKLNGTSAPV PGPRAQLKLN ASASDHRSNF SCSAALEIAG QVVHKHQGLE LHVLYGPRLD QRDCPGNWTW QEGSEQTLKC EAQGNPIPKL NCSRKGDGAS LPIGDLRPVR REVAGTYLCR ATSARGRVTR EVVLNVLHGQ NIL (453)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	50.0 kDa
<b>Entrez</b>	<a href="#">281839</a>
<b>Expiration Date</b>	6 months from date of receipt.