

Product Datasheet

Goat IL-1 alpha protein (orb1216039)

Description	The Caprine IL-1 alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	IL-1 alpha
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	SAHYSFQSNV KYNFMRVIHQ ECILNDALNQ SIIRDMSGPY LMAATLNNLE EAVKFDMVAY VSEEDSQLPV TLRISKQLF VSAQNEDEPV LLKEMPETPK IIKDETLLF FWEKHGSMY FKSVAHPKLF IATKQEKL VH MASGPPSITD FQILEK (156)
Purity	98%
Source	Yeast
MW	17.9 kDa
Entrez	100861289
Expiration Date	6 months from date of receipt.