

Product Datasheet

Human CCL11 protein (orb1216058)

Description	The Human CCL11 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	CCL11
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GPASVPTTCC FNLANRKIPL QRLESYRRIT SGKCPQKAVI FTKLAKDIC ADPKKKWVQD SMKYLDQKSP TPKP (74)
Purity	98%
Source	Yeast
MW	8.4 kDa
Entrez	6356
Expiration Date	6 months from date of receipt.