

Product Datasheet

Guinea pig LIF protein (orb1216074)

| | |
|-------------------------|---|
| Description | The Guinea Pig LIF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ... |
| Target | LIF |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | SPLPITPVNA TCATRHPCPS NVMSQIRNQL AQLNSSANAL FILYYTAQGE PFPNNLDKLC APNVTFDFPPF HANGTEKVKL VELYRIIAYL GAALGNITHD QKLLNPALS LHSKLSVTAD VMRGLLSNVH CRLCNKYHVA HVDVTYGPDT SDKDVFQKKK LGCQLLGKYK QVIGVVTQAF (180) |
| Purity | 98% |
| Source | Yeast |
| MW | 19.8 kDa |
| Entrez | 100713602 |
| Expiration Date | 6 months from date of receipt. |