

Product Datasheet

Goat IL-17A protein (orb1216096)

Description	The Caprine IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. ...
Target	IL-17A
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GVIIQSPGC PPTEDKNFPQ HVRVNLNIVN RNTNSRRPTN YHKRSTSPWT LHRNEDPERY PSVIWEAKCS HSGCINAEGK VDHHMNSVTI QQEILVLRRE SQHCPHSFRL EKMLVAVGCT CVTPIVRHVA (130)
Purity	98%
Source	Yeast
MW	14.8 kDa
Entrez	100860877
Expiration Date	6 months from date of receipt.