

## Product Datasheet

### Feline IL-17A protein (orb1216097)

<b>Description</b>	The Feline IL-17A yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T...
<b>Target</b>	IL-17A
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GIAFPQNPGC PTTEDKNFPQ HVKVNVLNLN GNKSSRRPLD YYRRSTSPWS LHRNEDPERY PSVIWEAKCL HWGCVNTEGK EDHHMNSVPI QQEILVLRRE SRHCPHSFRL EKMLVTVGCT CVTPIVRHVV (130)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	15.1 kDa
<b>Entrez</b>	<a href="#">101095339</a>
<b>Expiration Date</b>	6 months from date of receipt.