

## Product Datasheet

### Canine IFN beta 1 protein (orb1216100)

<b>Description</b>	The Canine IFN beta 1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IFN beta 1
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MSNDLLRSQL SSSSLECQEL LLQLNGTTEY CLKDRINFEI PEEIEKSRQF QKEDIILITH EMFQKIFDIF RRNISRTGWN ETTVENLLVK LHWQKEHLEI ILEDVKEKEN FTWDNRTLLH LKKYYLRIVQ YLKAKEYSIC AWTIVQAEIC RNFFFLNILT DYLQN (165)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	20.0 kDa
<b>Entrez</b>	<b>481558</b>
<b>Expiration Date</b>	6 months from date of receipt.