

## Product Datasheet

### Porcine IFN omega 5 protein (orb1216102)

<b>Description</b>	The Swine IFN omega 5 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
<b>Target</b>	IFN omega 5
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	GSLGCDLSQN HVHISRKNLV LLHQMRRISP SFCLKDRKDF GLPQEMVDGS HLQKAQAISV LHEMLQQTFL LFHIKRSSAA WDSILLDKLH SGLHQQLEDL DPCLVQEMGE QASALGMAMK KYFQGIHLYL KEKKYSDCAW EIVRVEIMRA LSFSTNLQER LRIMDEHLGS P (171)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	19.7 kDa
<b>Entrez</b>	<b>100154098</b>
<b>Expiration Date</b>	6 months from date of receipt.