

Product Datasheet

Monkey IFN gamma protein (orb1216104)

| | |
|-------------------------|--|
| Description | The Cynomolgus Monkey IFN gamma yeast-derived recombinant protein is not tagged and is naturally end... |
| Target | IFN gamma |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | QDPYVKEAEN LKKYFNAGDP DVADNGTLFL DILRNWKEES DRKIMQSQIV SFYFKLFKNF KDDQRIQKSV ETIKEDINVK FFNSNKKKRD DFEKLTNYSV TDSNVQRKAV HELIQVMAEL SPAAKIGKRK RSQMFRGRRA SQ (142) |
| Purity | 98% |
| Source | Yeast |
| MW | 16.8 kDa |
| Entrez | 102128291 |
| Expiration Date | 6 months from date of receipt. |