

Product Datasheet

Bovine CD40LG protein (orb1216113)

Description	The Bovine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MHKGDQEPQI AAHVISEASS KTTSVLQWAP KGYTLSNNL VTLENGKQLA VKRQGFYYIY TQVTFCSNRE TLSQAPFIAS LCLKSPSGSE RILLRAANTH SSSKPCGQQS IHLGGVFELQ SGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.1 kDa
Entrez	282387
Expiration Date	6 months from date of receipt.