

Product Datasheet

Human CD40LG protein (orb1216122)

Description	The Human & Monkey CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endo...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQKGDQNPQI AAHWISEASS KTTSVLQWAE KGYTMSNNL VTLENGKQLT VKRQGLYYIY AQVTFCSNRE ASSQAPFIAS LCLKSPGRFE RILLRAANTH SSAKPCGQQS IHLGGVFELQ PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.2 kDa
Entrez	959
Expiration Date	6 months from date of receipt.