

Product Datasheet

Feline CD40 Ligand protein (orb1216123)

Description	The Feline CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CD40 Ligand
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQKGDQDPRV AAHVISEASS STASVLQWAP KGYTISNL VTLENGKQLA VKRQGLYYIY AQVTFCSNRE ASSQAPFIAS LCLHSPSGSE RVLLRAANAR SSSKPCGQQS IHLGGVFELH PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	15.9 kDa
Entrez	493850
Expiration Date	6 months from date of receipt.