

Product Datasheet

Equine CD40LG protein (orb1216124)

Description	The Equine CD40 Ligand yeast-derived recombinant protein is not tagged and is naturally endotoxin-fr...
Target	CD40LG
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	MQKGDQEPQI AAHVISEASS KTASVLQWAQ KGYTISNNL VTLENGKQLA VKRQGLYYIY AQVTFCSNRE ASGQAPFIAS LCLRSVSGSE RILLRAANTH SSSKPCGQQS IHLGGVFELQ PGASVFNVT DPSQVSHGTG FTSFGLLKL (149)
Purity	98%
Source	Yeast
MW	16.0 kDa
Entrez	100054624
Expiration Date	6 months from date of receipt.