

Product Datasheet

Canine CX3CL1 protein (orb1216126)

| | |
|-------------------------|---|
| Description | The Canine CX3CL1 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. T... |
| Target | CX3CL1 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | QHLGVKKCNV TCHKMTSEIP VALLVHYQRN QESCGKPAII LKTKQNRIFC ADPKERWVQK AIAHLEHQIA VTIRNG (76) |
| Purity | 98% |
| Source | Yeast |
| MW | 8.7 kDa |
| Entrez | 487265 |
| Expiration Date | 6 months from date of receipt. |