

Product Datasheet

Rabbit IL-22 protein (orb1216135)

| | |
|-------------------------|---|
| Description | The Rabbit IL-22 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. Th... |
| Target | IL-22 |
| Form/Appearance | Lyophilized |
| Storage | -20°C |
| Note | For research use only |
| Protein Sequence | LPISSHCRLD ESNFTSYITN LTFMLAKEAN SADNNTDVRL IGNELFHGVN EKDH CYLLKH VLNFTLEEV LFPESDRFQPY MQQVLSFLAK LSNNLSQCHI EGDDQHIQRN VQQLKDTVKQ LGESGEIKAI GELDLLFMAL RTACVSPEQS WKMS (154) |
| Purity | 98% |
| Source | Yeast |
| MW | 17.6 kDa |
| Entrez | 100348980 |
| Expiration Date | 6 months from date of receipt. |