

Product Datasheet

Porcine IL-2 protein (orb1216146)

Description	The Swine IL-2 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The ...
Target	IL-2
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APTSSSTKNT KKQLEPLLLD LQLLKEVKN YENADLSRML TFKFYMPKQA TELKHLQCLV EELKALEGVL NLGQSKNSDS ANIKESMNNI NVTVLELKGS ETSFKCEYDD ETVTAVEFLN KWITFCQSIY STL (134)
Purity	98%
Source	Yeast
MW	15.2 kDa
Entrez	396868
Expiration Date	6 months from date of receipt.