

Product Datasheet

Porcine EPO protein (orb1216147)

Description	The Swine Erythropoietin yeast-derived recombinant protein is not tagged and is naturally endotoxin-...
Target	EPO
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	APPRLICDSR VLERYILEAK EGENATMGCA ESCSFSENI T VPDTKVNFYA WKRMEVQQQA MEVWQGLALL SEAILQGQAL LANSSQPSEA LQLHVDKAVS GLRSLTSLR ALGAQKEAIP LPDASPSSAT PLRTFAVDTL CKLFRNYSNF LRGKLTLYTG EACRRRDR (168)
Purity	98%
Source	Yeast
MW	18.6 kDa
Entrez	397249
Expiration Date	6 months from date of receipt.