

Product Datasheet

Porcine IL-21 protein (orb1216148)

Description	The Swine IL-21 yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	IL-21
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	HKSSFQGQDR LLIRLRQLID TVDQLKNYVH DLDPPELLPAP EDVQRHCEQS AFSCFQKVEL KSANTGDNEK IINVLIKQLK RKLPTNAGR RQKHGLTCPT CDSYEKKPIK EFLERLKSLI QKMIHQHLS (129)
Purity	98%
Source	Yeast
MW	15.0 kDa
Entrez	403123
Expiration Date	6 months from date of receipt.