

Product Datasheet

Rabbit BAFF protein (orb1216156)

Description	The Rabbit BAFF yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
Target	BAFF
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	AVEGVEETVI QDCLQLIADS DTPIIRKGSY TFVPWLLSFK RGRALEEKEN KIVVKETGYF FIYGQVLYTD STFAMGHLIQ RKKVHVFGDE LSLVTLFRCI QNMPETLPNN SCYSAGIAKL EEGDELQLAI PRENAQISRD GDGTFFGALK LL (152)
Purity	98%
Source	Yeast
MW	17.1 kDa
Entrez	100101577
Expiration Date	6 months from date of receipt.