

Product Datasheet

Equine IL-12 p40 protein (orb1216158)

Description	The Equine IL-12/IL-23 p40 yeast-derived recombinant protein is not tagged and is naturally endotoxi...
Target	IL-12 p40
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	IWELEKDVYV VELDWYPDAP GEMVVLTCNT PEEEGITWTS AQSNEVLGSG KTLTIQVKEF GDAGWYTCHK GGEVLSHSHL LLHKKEDGIW STDILKDQKE SKNKTFLKCE AKNYSGRFTC WWLTAISTDL KFSVKSSRGS SDPRGVTCGA ATLSAERVSV DDREYKKYTV ECLEGSACPA AEESLPIEIV VDAVHKLKYE NYTSGFFIRD IIKPDPPKNL QLKPLKNSRQ VEVSWEYPET WSTPHSYFSL TFSIQVQGKN KKERKDRLFM DETSATVTCH KDGQIRVQAR DRYSSSWSE WASVSCS (307)
Purity	98%
Source	Yeast
MW	34.9 kDa
Entrez	100034222
Expiration Date	6 months from date of receipt.