

## Product Datasheet

### Sheep TNF alpha protein (orb1216169)

<b>Description</b>	The Ovine TNF alpha yeast-derived recombinant protein is not tagged and is naturally endotoxin-free....
<b>Target</b>	TNF alpha
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	LRSSSQASNN KPVAHVANI SAPGQLRWGD SYANALMANG VELKDNQLVV PTDGLYLIYS QVLFRGHGCP STPLFLTHTI SRIAVSYQTK VNILSAIKSP CHRETLEGAE AKPWYEPIYQ GGVFQLEKGD RLSAEINLPE YLDYAESGQV YFGIIAL (157)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.2 kDa
<b>Entrez</b>	<a href="#">443540</a>
<b>Expiration Date</b>	6 months from date of receipt.