

## Product Datasheet

### Mouse TWEAK protein (orb1216173)

<b>Description</b>	The Mouse TWEAK yeast-derived recombinant protein is not tagged and is naturally endotoxin-free. The...
<b>Target</b>	TWEAK
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Protein Sequence</b>	SAPKGRKARP RRAIAAHYEV HPRPGDGAQ AGVDGTVSGW EETKINSSSP LRYDRQIGEF TVIRAGLYYL YCQVHFDEGK AVYKLDLLV NGVLALRCLE EFSATAASSP GPQLRLCQVS GLLPLRPGSS LRIRTLPSWAH LKAAPFLTYF GLFQVH (156)
<b>Purity</b>	98%
<b>Source</b>	Yeast
<b>MW</b>	17.2 kDa
<b>Entrez</b>	<b>21944</b>
<b>Expiration Date</b>	6 months from date of receipt.